


MedicCleanAir®

- Swiss Made -

Information for the Patient



Our vision is an infection free Hospital

We see

- protected patients and hospital staff
- professional isolation (A.I.I.)
- rapid deployment at low cost
- clean air + negative pressure
- protection against Tuberculosis, MRSA, SARS etc.

We see

- professional isolation (P.E.)
- rapid deployment at low cost
- clean air + positive pressure
- protection for vulnerable patients like transplantation, neonates, burn units etc.

We see

- protection during building activities
- rapid deployment at low cost
- during small and large activities
- protection against aspergillus outbreaks, protection against epidemics

Proven Outcomes - Is what we deliver right now. **Outcomes** that result from clinical studies. **Outcomes** that dramatically improve your efforts in infection control. **Outcomes** that lead to a level of care that is exceptional.

Proof positive
„break the chain“

What is your vision?

TheHospital not only wants to give the patients the best medical service possible, but also the guarantee of an optimal protection through preventive hygiene measures.

Clean Hepa Filtered Air is one of these hygiene measures that have been implemented in your department.

TheHospital has decided to install MedicCleanAir-Pro Mobile Hepa Units in this department. These MedicCleanAir-Pro Units already function very successfully not only in various countries in Europe like Germany, United Kingdom, Italy, Belgium, and in Scandinavia, but in the Middle East and Africa as well.

The units are the only mobile air cleaning systems that are capable to effectively reduce airborne contamination, more important: reduce airborne infections by providing an optimal air quality.

The MedicCleanAir Concept has been thoroughly tested and evaluated by independent institutions, laboratories, independent University Hospitals and others – yours included -. These medical studies have been world-widely accepted and published in medical magazines by esteemed organisations.

The units clean the air by sucking the air of the patient's room into the unit over 360°, purify the air and return it back to the room over 360°. The fastest way to clean indoor air.

The filter cartridge is manufactured out of 3 filter systems (F9+activated carbon+Hepa14). This cartridge assures an optimal air-quality in the room.

The filter cartridge adsorbs particles, gases/odours and a combination of both.

Your help and support for the implementation of these infection control units is of the highest importance, not only for yourself, but also for the medical staff, visitors and other patients.

THEREFORE – ATTENTION !!

**The units must function continuously (24/24 hours) at the spot where they are installed.
Only medical staff is allowed to implement any changes to the operation of the units.**

We thank you on beforehand for your support to implement this new hospital hygiene measure that is very important - not in the least for you personally -.

In case you have any questions: the medical staff is at your disposal.

In case you have any technical questions regarding the MedicCleanAir Mobile Hepa Unit, do not hesitate to contact us.

Tel: +32-38607060 – Email: info@mediccleanair.com – Url: www.mediccleanair.com